



Product specification

Filtropur V50, Vacuum filtration unit, 500 ml, PES, 0.1 µm



Product information

Product number	83.3941.502
----------------	-------------

Product descriptions

Name	Filtropur V50, Vacuum filtration unit, 500 ml, PES, 0.1 µm
Product number	83.3941.502
Product description	Filtropur V50, Vacuum filtration unit, 500 ml, screw cap, membrane: PES, Ø membraner: 75 mm, pore size: 0.1 µm, for removing mycoplasmas, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

Material & colours

Color of product	transparent
Membrane material	Polyethersulfone (PES)
Closure material	High Density Polyethylene (HD-PE)

Size

Height of product	24.6 cm
Sample volume	500 ml
Filtration volume	500 ml

Product characteristics

Filter	yes
--------	-----

Area of the filter	38.04 cm ²
Graduation	yes
Pore size	0.1 µm
Closure type	Screw cap
Connection	Step taper
Example of use	for removing mycoplasmas

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / case	12
Piece(s) / pallet	288
Height of case	285 mm
Width of case	391 mm
Depth of case	508 mm
Volume of case	0.0566 cbm
Weight of case	3.8 kg
EAN inner box	4038917969628
EAN case	4038917969635

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT