




Product specification

Pipette tip, 1,250 µl, transparent, Biosphere® plus, 96 piece(s)/box



Product information

Product number	70.3070.205
----------------	-------------

Product descriptions

Name	Pipette tip, 1,250 µl, transparent, Biosphere® plus, 96 piece(s)/box
Product number	70.3070.205
Product description	Pipette tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere® plus, 96 piece(s)/box

Material & colours

Color of product	transparent
------------------	-------------

Size

Diameter	9 mm
Length of product	108 mm
Max. volume	1250 µl

Product characteristics

Filter	no
Graduation	yes
Autoclavable (max. °C)	121 °C

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	46080
Height of inner box	146 mm
Width of inner box	140 mm
Depth of inner box	520 mm
Height of case	315 mm
Width of case	298 mm
Depth of case	539 mm
Volume of case	0.0506 cbm
Weight of case	5.403 kg
EAN inner box	4066916022914
EAN case	4066916022921

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT