




## Product specification

**Pipette tip, 1,250 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill**



## Product information

|                |             |
|----------------|-------------|
| Product number | 70.3070.305 |
|----------------|-------------|

## Product descriptions

|                     |  |
|---------------------|--|
| Name                | Pipette tip, 1,250 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill                                      |
| Product number      | 70.3070.305  |
| Product description | Pipette tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere® plus, 96 piece(s)/SingleRefill |

## Material & colours

|                  |             |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

## Size

|                   |         |
|-------------------|---------|
| Diameter          | 9 mm    |
| Length of product | 108 mm  |
| Max. volume       | 1250 µl |

## Product characteristics

|            |     |
|------------|-----|
| Filter     | no  |
| Graduation | yes |

|                        |        |
|------------------------|--------|
| Autoclavable (max. °C) | 121 °C |
|------------------------|--------|

## Purity & certification

|                  |                                      |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD)            |
| CE certified     | CE - manufacturer's self-declaration |
| Purity standard  | Biosphere plus                       |
| Sterilisation    | Gas sterilisation                    |
| Batched          | yes                                  |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 480           |
| Piece(s) / case      | 1920          |
| Piece(s) / pallet    | 57600         |
| Height of inner box  | 146 mm        |
| Width of inner box   | 130 mm        |
| Depth of inner box   | 475 mm        |
| Height of case       | 315 mm        |
| Width of case        | 278 mm        |
| Depth of case        | 494 mm        |
| Volume of case       | 0.0433 cbm    |
| Weight of case       | 4.53 kg       |
| EAN inner box        | 4066916024352 |
| EAN case             | 4066916024369 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT