




## Product specification

**Filter tip, 1,250 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill**



## Product information

Product number	70.3070.355
----------------	-------------

## Product descriptions

Name	Filter tip, 1,250 µl, transparent, Biosphere® plus, 96 piece(s)/SingleRefill
Product number	70.3070.355
Product description	Filter tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere® plus, 96 piece(s)/SingleRefill

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	9 mm
Length of product	108 mm
Max. volume	1250 µl

## Product characteristics

Filter	yes
Graduation	yes

Autoclavable (max. °C)	121 °C
------------------------	--------

## Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere plus
Sterilisation	Gas sterilisation
Batched	yes

## Packaging

Piece(s) / inner box	480
Piece(s) / case	1920
Piece(s) / pallet	57600
Height of inner box	146 mm
Width of inner box	130 mm
Depth of inner box	475 mm
Height of case	315 mm
Width of case	278 mm
Depth of case	494 mm
Volume of case	0.0433 cbm
Weight of case	4.53 kg
EAN inner box	4066916025144
EAN case	4066916025151

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT