



## Product specification

**Sarpette® M, single-channel, variable, 0.5 µl - 10 µl, grey, 1 piece(s)/case**



### Product information

Product number	90.3100.012
----------------	-------------

### Product descriptions

Name	Sarpette® M, single-channel, variable, 0.5 µl - 10 µl, grey, 1 piece(s)/case
Product number	90.3100.012
Product description	Sarpette® M, single-channel, variable, 0.5 µl - 10 µl, grey, 1 piece(s)/case

### Material & colours

Color of product	grey
------------------	------

### Product characteristics

Autoclavable (max. °C)	121 °C
------------------------	--------

### Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	no

### Packaging

Height of inner box	30 mm
---------------------	-------

Width of inner box	90 mm
Depth of inner box	315 mm
Height of case	30 mm
Width of case	90 mm
Depth of case	320 mm
Volume of case	0.0009 cbm
Weight of case	0.18 kg
EAN inner box	4066916002633
EAN case	4066916002640

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT