



## Product specification

**Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP, with paper label**



## Product information

Product number	63.595.254
----------------	------------

## Product descriptions

Name	Screw cap tube, 25 ml, (LxØ): 90 x 25 mm, PP, with paper label
Product number	63.595.254
Product description	Screw cap tube, working volume: 25 ml, (LxØ): 90 x 25 mm, material: PP, skirted conical base, transparent, screw cap, yellow, cap assembled, with paper label, label/print: white/yellow, Patient information, with graduations, sterile, 500 piece(s)/bag

## Material & colours

Color of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Color of cap	yellow

## Size

Diameter	25 mm
Length of product	90 mm

## Product characteristics

Base shape	skirted conical base
Label/Print	With paper label
Color of print/label	white/yellow
Graduation	yes
Closure type	Screw cap
Centrifugation max (RCF)	4500 x g

## Purity & certification

Satisfies the requirement	IATA
Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

## Packaging

Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	15000
Height of case	260 mm
Width of case	408 mm
Depth of case	458 mm
Volume of case	0.0486 cbm
Weight of case	5.0466 kg
EAN inner box	4066916009717
EAN case	4066916009724

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT