



Product specification

Filtropur BT 25, Bottle top filter, 250 ml, PES, 0.22 µm



Product information

Product number	83.3940.511
----------------	-------------

Product descriptions

Name	Filtropur BT 25, Bottle top filter, 250 ml, PES, 0.22 µm
Product number	83.3940.511
Product description	Filtropur BT 25, Bottle top filter, for screw cap 45 mm, 250 ml, membrane: PES, Ø membraner: 50 mm, pore size: 0.22 µm, for sterile filtration, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

Material & colours

Color of product	transparent
Membrane material	Polyethersulfone (PES)

Size

Height of product	11.35 cm
Sample volume	250 ml
Filtration volume	250 ml

Product characteristics

Filter	yes
Area of the filter	15.16 cm²
Graduation	yes

Pore size	0.22 µm
Connection	Step taper
Example of use	for sterile filtration

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / case	24
Piece(s) / pallet	840
Height of case	236 mm
Width of case	358 mm
Depth of case	433 mm
Volume of case	0.0366 cbm
Weight of case	2.8 kg
EAN inner box	4038917969406
EAN case	4038917969512

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

