



Product specification

Collection vessel, 1,000 ml



Product information

Product number	83.3942.505
----------------	-------------

Product descriptions

Name	Collection vessel, 1,000 ml
Product number	83.3942.505
Product description	Collection vessel, for Filtropur BT, 1,000 ml, screw cap, cap assembled, sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase, 1 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	Polypropylene (PP)

Size

Height of product	17.57 cm
Sample volume	1000 ml
Filtration volume	1000 ml

Product characteristics

Filter	no
Vessel type	Flask
Graduation	yes
Closure type	Screw cap

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic, non-cytotoxic, free from DNase/Rnase
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / case	24
Piece(s) / pallet	384
Height of case	385 mm
Width of case	425 mm
Depth of case	545 mm
Volume of case	0.0892 cbm
Weight of case	6.09 kg
EAN inner box	4038917969178
EAN case	4038917969185

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

