


# Product specification

**Filter tip, 1,250 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/box**



## Product information

|                |             |
|----------------|-------------|
| Product number | 70.3070.275 |
|----------------|-------------|

## Product descriptions

|      |  |
|------|--|
| Name | Filter tip, 1,250 µl, transparent, Biosphere® plus, Low retention, 96 piece(s)/box |
|------|--|

|                |             |
|----------------|-------------|
| Product number | 70.3070.275 |
|----------------|-------------|

|                     |   |
|---------------------|---|
| Product description | Filter tip, working volume: 1,250 µl, transparent, filling level rings, Biosphere® plus, Low retention, 96 piece(s)/box |
|---------------------|---|

## Material & colours

|                  |             |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

## Size

|          |      |
|----------|------|
| Diameter | 9 mm |
|----------|------|

|                   |        |
|-------------------|--------|
| Length of product | 108 mm |
|-------------------|--------|

|             |         |
|-------------|---------|
| Max. volume | 1250 µl |
|-------------|---------|

## Product characteristics

|        |     |
|--------|-----|
| Filter | yes |
|--------|-----|

|            |     |
|------------|-----|
| Graduation | yes |
|------------|-----|

|                        |        |
|------------------------|--------|
| Autoclavable (max. °C) | 121 °C |
|------------------------|--------|

## Purity & certification

|                  |                                      |
|------------------|--------------------------------------|
| Product category | In-vitro diagnostic (IVD)            |
| CE certified     | CE - manufacturer's self-declaration |
| Purity standard  | Biosphere plus                       |
| Sterilisation    | Gas sterilisation                    |
| Batched          | yes                                  |

## Packaging

|                      |               |
|----------------------|---------------|
| Piece(s) / inner box | 480           |
| Piece(s) / case      | 1920          |
| Piece(s) / pallet    | 46080         |
| Height of inner box  | 146 mm        |
| Width of inner box   | 140 mm        |
| Depth of inner box   | 520 mm        |
| Height of case       | 315 mm        |
| Width of case        | 298 mm        |
| Depth of case        | 539 mm        |
| Volume of case       | 0.0506 cbm    |
| Weight of case       | 5.403 kg      |
| EAN inner box        | 4066916031992 |
| EAN case             | 4066916032005 |

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

