

Product specification

Cell culture plate, 6 well, surface: BIOFLOAT™, flat base



Product information

| | |
|----------------|-------------|
| Product number | 83.3920.400 |
|----------------|-------------|

Product descriptions

| | |
|------|---|
| Name | Cell culture plate, 6 well, surface: BIOFLOAT™, flat base |
|------|---|

| | |
|----------------|-------------|
| Product number | 83.3920.400 |
|----------------|-------------|

| | |
|---------------------|---|
| Product description | Cell culture plate, 6 well, material: PS, surface: BIOFLOAT™, for spheroid culture, flat base, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag |
|---------------------|---|

Material & colours

| | |
|------------------|-------------|
| Color of product | transparent |
|------------------|-------------|

| | |
|------------------|------------------|
| Closure material | Polystyrene (PS) |
|------------------|------------------|

Size

| | |
|------------------|---------|
| Width of product | 85.2 mm |
|------------------|---------|

| | |
|-------------------|----------|
| Length of product | 127.8 mm |
|-------------------|----------|

| | |
|-------------------|---------|
| Height of product | 22.4 mm |
|-------------------|---------|

| | |
|-------------|----------|
| Max. volume | 15.75 ml |
|-------------|----------|

Product characteristics

| | |
|-------------|----------------------|
| Growth area | 8.87 cm ² |
|-------------|----------------------|

| | |
|------------|-----------|
| Base shape | Flat base |
|------------|-----------|

| | |
|--------------|------------------|
| Vessel type | Plate |
| Closure type | Slip-on lid |
| Cell type | spheroid culture |

Purity & certification

| | |
|------------------|--|
| Product category | not regulated |
| CE certified | no CE mark |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Electron irradiation |
| Batched | yes |

Packaging

| | |
|---------------------|---------------|
| Piece(s) / case | 18 |
| Piece(s) / pallet | 1080 |
| Height of inner box | 80 mm |
| Width of inner box | 140 mm |
| Depth of inner box | 205 mm |
| Height of case | 265 mm |
| Width of case | 257 mm |
| Depth of case | 322 mm |
| Volume of case | 0.0219 cbm |
| Weight of case | 1.9225 kg |
| EAN inner box | 4066916042783 |
| EAN case | 4066916042790 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.

