




Product specification

Mailing envelope, PE, (LxW): 310 x 123 mm



Product information

| | |
|----------------|---------|
| Product number | 95.1411 |
|----------------|---------|

Product descriptions

| | |
|---------------------|-----------------------------------------------------------------------------|
| Name | Mailing envelope, PE, (LxW): 310 x 123 mm |
| Product number | 95.1411 |
| Product description | Mailing envelope, material: PE, blue, (LxW): 310 x 123 mm, 50 piece(s)/case |

Material & colours

| | |
|------------------|-------------------|
| Product material | Polyethylene (PE) |
| Color of product | blue |

Size

| | |
|-------------------|--------|
| Width of product | 123 mm |
| Length of product | 310 mm |

Product characteristics

| | |
|----------------|----|
| Suction insert | no |
|----------------|----|

Purity & certification

| | |
|---------------------------|-------------------------|
| Satisfies the requirement | P 650 |
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|------------------------|---------------|
| Minimum order quantity | 1000 |
| Piece(s) / inner box | 50 |
| Piece(s) / case | 1000 |
| Piece(s) / pallet | 40000 |
| Height of case | 130 mm |
| Width of case | 323 mm |
| Depth of case | 428 mm |
| Volume of case | 0.018 cbm |
| Weight of case | 8.42 kg |
| EAN inner box | 4038917086547 |
| EAN case | 4038917005098 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT