




## Product specification

**Mailing envelope, PE, (LxW): 310 x 123 mm**



### Product information

Product number	95.1411
----------------	---------

### Product descriptions

Name	Mailing envelope, PE, (LxW): 310 x 123 mm
Product number	95.1411
Product description	Mailing envelope, material: PE, blue, (LxW): 310 x 123 mm, 50 piece(s)/case

### Material & colours

Color of product	blue
------------------	------

### Size

Width of product	123 mm
Length of product	310 mm

### Product characteristics

Suction insert	no
----------------	----

### Purity & certification

Satisfies the requirement	P 650
Product category	not regulated
CE certified	no CE mark
Batched	no

## Packaging

Piece(s) / inner box	50
Piece(s) / case	1000
Piece(s) / pallet	40000
Height of inner box	20 mm
Width of inner box	130 mm
Depth of inner box	317 mm
Height of case	173 mm
Width of case	355 mm
Depth of case	390 mm
Volume of case	0.024 cbm
Weight of case	8.42 kg
EAN inner box	4038917086547
EAN case	4038917005098

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT