



## Product specification

**flexiPERM®**, con A, re-usable silicon insert, 1 well, 5 piece(s)/bag



## Product information

Product number	94.6077.434
----------------	-------------

## Product descriptions

Name	flexiPERM®, con A, re-usable silicon insert, 1 Well, 5 piece(s)/Bag
Product number	94.6077.434
Product description	flexiPERM®, con A, re-usable silicon insert, for dividing the growth surface, 1 well, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

## Material & colours

Product material	silicone
Color of product	transparent

## Size

Diameter	12 mm
Height of product	8 mm

## Product characteristics

Growth area	1.1 cm <sup>2</sup>
-------------	---------------------

Product category	no medical device   IVD
Purity standard	non-pyrogenic/endotoxin-free, non-cytotoxic
Batched	yes

## Packaging

Piece(s) / case	100
Piece(s) / pallet	7200
Height of case	225 mm
Width of case	253 mm
Depth of case	388 mm
Volume of case	0.0221 cbm
Weight of case	1.42 kg
EAN inner box	4038917255264
EAN case	4038917255257

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT