



## Product specification

**flexiPERM®, micro 12, re-usable silicon insert, 12 well, 5 piece(s)/bag**



## Product information

Product number	94.6011.436
----------------	-------------

## Product descriptions

Name	flexiPERM®, micro 12, re-usable silicon insert, 12 Well, 5 piece(s)/Bag
Product number	94.6011.436
Product description	flexiPERM®, micro 12, re-usable silicon insert, for dividing the growth surface, 12 well, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Width of product	25 mm
Length of product	68 mm
Height of product	10 mm

## Product characteristics

Growth area	0.3 cm <sup>2</sup>
Number of wells	12

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	non-pyrogenic/endotoxin-free, non-cytotoxic
Batched	yes

## Packaging

Piece(s) / case	100
Piece(s) / pallet	7200
Height of inner box	30 mm
Width of inner box	100 mm
Depth of inner box	185 mm
Height of case	225 mm
Width of case	253 mm
Depth of case	388 mm
Volume of case	0.0221 cbm
Weight of case	2.2 kg
EAN inner box	4038917254489
EAN case	4038917254472

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT