



## Product specification

**Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm and 12 x 75 mm tubes**



## Product information

Product number	83.3945.999
----------------	-------------

## Product descriptions

Name	Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm and 12 x 75 mm tubes
Product number	83.3945.999
Product description	Cell strainer adapter, for 15 ml, 5 ml, 17 x 100 mm and 12 x 75 mm tubes, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

## Material & colours

Product material	Polypropylene (PP)
Color of product	transparent

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

## Packaging

Minimum order quantity	25
Piece(s) / inner box	25
Piece(s) / case	200
Piece(s) / pallet	8400
Height of case	584 mm
Width of case	330 mm
Depth of case	229 mm
Volume of case	0.0441 cbm
Weight of case	2 kg
EAN inner box	4038917415323
EAN case	4038917415316

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT