




Product specification

Cell strainer, pore size: 100 µm, yellow



Product information

Product number	83.3945.100
----------------	-------------

Product descriptions

Name	Cell strainer, Pore size: 100 µm, yellow
Product number	83.3945.100
Product description	Cell strainer, pore size: 100 µm, yellow, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Filter material	Nylon
Membrane material	Nylon
Frame material	Polypropylene (PP)
Color of frame	yellow

Size

Diameter	34.09 mm
Height of product	21 mm

Product characteristics

Filter	yes
Area of the filter	4.2 cm ²
Pore size	100 µm

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Electron irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	50
Piece(s) / inner box	50
Piece(s) / case	600
Piece(s) / pallet	9600
Height of case	410 mm
Width of case	385 mm
Depth of case	520 mm
Volume of case	0.0821 cbm
Weight of case	6.78 kg
EAN inner box	4038917415040
EAN case	4038917415033

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT