



Product specification

Cell culture tubes, (LxØ): 125 x 16 mm, round base, TC-treated



Product information

Product number	83.9923.945
----------------	-------------

Product descriptions

Name	Cell culture tubes, (LxØ): 125 x 16 mm, Round base, TC-treated
Product number	83.9923.945
Product description	Cell culture tubes, 15 ml, (LxØ): 125 x 16 mm, round base, PS, TC-treated, for Suspension and mono-layer cultures, sterile, 5 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Closure material	Polyethylene (PE)
Color of cap	red

Size

Diameter	16 mm
Length of product	125 mm

Product characteristics

Base shape	Round base
------------	------------

Vessel type	Tube
Graduation	no
Type of surface treatment	TC-treated
Cap	standard screw cap
Cell type	Suspension and mono-layer cultures

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Minimum order quantity	1000
Piece(s) / case	1000
Piece(s) / pallet	21000
Height of case	372 mm
Width of case	438 mm
Depth of case	428 mm
Volume of case	0.0697 cbm
Weight of case	9.4 kg
EAN inner box	4038917154635
EAN case	4038917027892

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT