



Product specification

Cell culture flask, T-175, surface: Suspension, Filter cap



Product information

Product number	83.3912.502
----------------	-------------

Product descriptions

Name	Cell culture flask, T-175, Surface: Suspension, Filter cap
Product number	83.3912.502
Product description	Cell culture flask, T-175, material: PS, surface: Suspension, for suspension cells, coding colour: green, Filter cap, TC Tested, 5 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Closure material	Polyethylene (PE)
Color of cap	green

Product characteristics

Growth area	175 cm ²
Vessel type	Flask
Cap	Filter cap
Cell type	suspension cells

Purity & certification

Product category	no medical device IVD
------------------	-------------------------

Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	40
Piece(s) / case	40
Piece(s) / pallet	960
Height of case	450 mm
Width of case	267 mm
Depth of case	522 mm
Volume of case	0.0627 cbm
Weight of case	5.62 kg
EAN inner box	4038917337922
EAN case	4038917337915

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT