



Product specification

Cell culture flask, T-25, surface: Suspension, filter cap



Product information

Product number	83.3910.502
----------------	-------------

Product descriptions

Name	Cell culture flask, T-25, Surface: Suspension, Filter cap
Product number	83.3910.502
Product description	Cell culture flask, T-25, material: PS, surface: Suspension, for suspension cells, coding colour: green, Filter cap, TC Tested, 10 piece(s)/bag

Material & colours

Closure material	Polyethylene (PE)
Color of cap	green

Product characteristics

Growth area	24.7 cm ²
Vessel type	Flask
Cap	Filter cap
Cell type	suspension cells

Purity & certification

Product category	not regulated
CE certified	no CE mark

Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	10
Piece(s) / case	300
Piece(s) / pallet	6000
Height of case	350 mm
Width of case	408 mm
Depth of case	595 mm
Volume of case	0.085 cbm
Weight of case	7.5099 kg
EAN inner box	4038917337809
EAN case	4038917337793

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT