




Product specification

Tissue culture dish, (ØxH): 100 x 20 mm, surface:
Suspension



Product information

Product number	83.3902.500
----------------	-------------

Product descriptions

Name	Tissue culture dish, (ØxH): 100 x 20 mm, Surface: Suspension
Product number	83.3902.500
Product description	Tissue culture dish, (ØxH): 100 x 20 mm, material: PS, surface: Suspension, for suspension cells, coding colour: green, TC Tested, 10 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Closure material	Polystyrene (PS)

Size

Diameter	100 mm
Height of product	20 mm

Product characteristics

Growth area	58.8 cm ²
Base shape	SUREGrip Design

Vessel type	Dish
Closure type	Slip-on lid
Cell type	suspension cells

Purity & certification

Product category	no medical device IVD
Purity standard	TC Tested
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	300
Piece(s) / inner box	10
Piece(s) / case	300
Piece(s) / pallet	4200
Height of case	460 mm
Width of case	313 mm
Depth of case	505 mm
Volume of case	0.0727 cbm
Weight of case	7.3701 kg
EAN inner box	4038917338226
EAN case	4038917338219

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT