




Product specification

Tissue culture dish, (ØxH): 60 x 15 mm, surface:
Suspension



Product information

Product number	83.3901.500
----------------	-------------

Product descriptions

Name	Tissue culture dish, (ØxH): 60 x 15 mm, Surface: Suspension
Product number	83.3901.500
Product description	Tissue culture dish, (ØxH): 60 x 15 mm, material: PS, surface: Suspension, for suspension cells, coding colour: green, TC Tested, 10 piece(s)/bag

Material & colours

Closure material	Polystyrene (PS)
------------------	------------------

Size

Diameter	60 mm
Height of product	15 mm

Product characteristics

Growth area	22.1 cm ²
Base shape	SUREGrip Design
Vessel type	Dish

Closure type

Slip-on lid

Cell type

suspension cells

Purity & certification

Product category

no medical device | IVD

Purity standard

TC Tested

Sterilisation

Gamma irradiation

Pyrogen-free

yes

Batched

yes

Packaging

Piece(s) / inner box

10

Piece(s) / case

500

Piece(s) / pallet

12000

Height of case

380 mm

Width of case

332 mm

Depth of case

332 mm

Volume of case

0.0419 cbm

Weight of case

5.02 kg

EAN inner box

4038917338134

EAN case

4038917338127

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT