




## Product specification

**Tissue culture dish, (ØxH): 35 x 10 mm, surface:**  
**Suspension**



## Product information

Product number	83.3900.500
----------------	-------------

## Product descriptions

Name	Tissue culture dish, (ØxH): 35 x 10 mm, Surface: Suspension
Product number	83.3900.500
Product description	Tissue culture dish, (ØxH): 35 x 10 mm, material: PS, surface: Suspension, for suspension cells, coding colour: green, TC Tested, 10 piece(s)/bag

## Material & colours

Closure material	Polystyrene (PS)
------------------	------------------

## Size

Diameter	35 mm
Height of product	10 mm

## Product characteristics

Growth area	9.4 cm <sup>2</sup>
Base shape	SUREGrip Design
Vessel type	Dish

Closure type

Slip-on lid

Cell type

suspension cells

## Purity & certification

Product category

not regulated

CE certified

no CE mark

Purity standard

TC Tested

Sterilisation

Gamma irradiation

Pyrogen-free

yes

Batched

yes

## Packaging

Piece(s) / inner box

10

Piece(s) / case

500

Piece(s) / pallet

30000

Height of case

255 mm

Width of case

257 mm

Depth of case

322 mm

Volume of case

0.0211 cbm

Weight of case

2.59 kg

EAN inner box

4038917338011

EAN case

4038917365819

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT