



Product specification

Coverslips, Ø: 13 mm



Product information

Product number	83.1840.002
----------------	-------------

Product descriptions

Name	Coverslips, Ø: 13 mm
Product number	83.1840.002
Product description	Coverslips, for adherent cells, Ø: 13 mm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 25 piece(s)/bag

Material & colours

Product material	Polyethylene terephthalate, Glycol-modified (PET-G)
Color of product	transparent

Size

Diameter	13 mm
----------	-------

Product characteristics

Type of surface treatment	TC-treated
Cell type	adherent cells

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic

Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	200
Piece(s) / inner box	200
Piece(s) / case	2000
Piece(s) / pallet	288000
Height of case	335 mm
Width of case	144 mm
Depth of case	310 mm
Volume of case	0.015 cbm
Weight of case	1.8 kg
EAN inner box	4038917134057
EAN case	4038917025393

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT