



Product specification

Filtropur BT 100, Bottle top filter, 1,000 ml, PES, 0.2 µm



Product information

Product number	83.3942.101
----------------	-------------

Product descriptions

Name	Filtropur BT 100, Bottle top filter, 1,000 ml, PES, 0.2 µm
Product number	83.3942.101
Product description	Filtropur BT 100, Bottle top filter, for screw cap 45 mm, 1,000 ml, membrane: PES, Ø membraner: 90 mm, pore size: 0.2 µm, for sterile filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Material & colours

Color of product	transparent
Membrane material	Polyethersulfone (PES)

Size

Height of product	15.27 cm
Sample volume	1000 ml
Filtration volume	1000 ml

Product characteristics

Filter	yes
Area of the filter	63.6 cm ²

Graduation	yes
Pore size	0.2 µm
Connection	Step taper
Example of use	for sterile filtration

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	24
Piece(s) / pallet	288
Height of case	439 mm
Width of case	343 mm
Depth of case	572 mm
Volume of case	0.0861 cbm
Weight of case	5.08 kg
EAN inner box	4038917412353
EAN case	4038917412346

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT