



Product specification

Collection vessel, 500 ml



Product information

Product number	83.3941.005
----------------	-------------

Product descriptions

Name	Collection vessel, 500 ml
Product number	83.3941.005
Product description	Collection vessel, for Filtropur BT 50, 500 ml, screw cap, cap assembled, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Closure material	Polypropylene (PP)

Size

Height of product	13.21 cm
Sample volume	500 ml
Filtration volume	500 ml

Product characteristics

Filter	no
Vessel type	Flask
Graduation	yes
Closure type	Screw cap

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	24
Piece(s) / case	24
Piece(s) / pallet	288
Height of case	444 mm
Width of case	343 mm
Depth of case	572 mm
Volume of case	0.0871 cbm
Weight of case	3.63 kg
EAN inner box	4038917412414
EAN case	4038917412407

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT