

# Product specification

## Collection vessel, 500 ml



### Product information

Order number	83.3941.005
--------------	-------------

### Product descriptions

Name	Collection vessel, 500 ml
Product number	83.3941.005
Product description	Collection vessel, for Filtropur BT 50, 500 ml, screw cap, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

### Product characteristics

Filter	no
Tube	Flask
Graduation	yes
Type of cap	Screw cap
Suitable for	Filtropur BT 50

### Size

Height of product	13.21 cm
Sample volume	500 ml
Filtration volume	500 ml

### Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Cap material	Polypropylene (PP)

### Purity & certification

Accordance with MDD or IVDD	no medical device   IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic

Sterilisation	Gamma irradiation
Pyrogen-free	yes
Lot or serial number	yes

## Packaging

Piece(s) / case	24
Piece(s) / pallet	288
Volume of case	0.0873 cbm
Weight of case	3.63 kg
EAN inner box	4038917412414
EAN case	4038917412407

**This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.**

**This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.**

This document was prepared by EDP support and is valid without signature.

