




Product specification

Collection vessel, 250 ml



Product information

Product number	83.3940.005
----------------	-------------

Product descriptions

Name	Collection vessel, 250 ml
Product number	83.3940.005
Product description	Collection vessel, for Filtropur BT 25, 250 ml, screw cap, cap assembled, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Material & colours

Color of product	transparent
Closure material	Polypropylene (PP)

Size

Diameter	92.1 mm
Height of product	88.3 mm
Sample volume	250 ml
Filtration volume	250 ml

Product characteristics

Filter	no
Vessel type	Flask
Graduation	yes
Closure type	Screw cap

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	24
Piece(s) / pallet	576
Height of case	381 mm
Width of case	267 mm
Depth of case	521 mm
Volume of case	0.0529 cbm
Weight of case	2.72 kg
EAN inner box	4038917412384
EAN case	4038917412377

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT