



Product specification

Syringe filter, Filtropur S plus, CA, pore size: 0.2 µm, for sterile filtration



Product information

Product number	83.1826.102
----------------	-------------

Product descriptions

Name	Syringe filter, Filtropur S plus, CA, Pore size: 0.2 µm, for sterile filtration
Product number	83.1826.102
Product description	Syringe filter, Filtropur S plus, pre-filter: GF, membrane: CA, filtration surface: 6.2 cm², pore size: 0.2 µm, for sterile filtration, connection: Luer lock, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister

Material & colours

Color of product	light blue
Membrane material	Cellulose acetate (CA)

Size

Diameter	33 mm
----------	-------

Product characteristics

Filter	yes
Area of the filter	6.2 cm²
Pore size	0.2 µm

Connection	Luer lock
Example of use	for sterile filtration

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	12000
Height of inner box	107 mm
Width of inner box	133 mm
Depth of inner box	200 mm
Height of case	280 mm
Width of case	365 mm
Depth of case	525 mm
Volume of case	0.0537 cbm
Weight of case	4.2 kg
EAN inner box	4038917133920
EAN case	4038917025263

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT