



Product specification

Filtropur V100, Vacuum filtration unit, 1,000 ml, PES, 0.22 µm



Product information

Product number	83.3942.001
----------------	-------------

Product descriptions

Name	Filtropur V100, Vacuum filtration unit, 1,000 ml, PES, 0.22 µm
Product number	83.3942.001
Product description	Filtropur V100, Vacuum filtration unit, 1,000 ml, screw cap, membrane: PES, Ø membraner: 90 mm, pore size: 0.22 µm, for sterile filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent
Membrane material	Polyethersulfone (PES)
Closure material	Polypropylene (PP)

Size

Height of product	28.33 cm
Sample volume	1000 ml
Filtration volume	1000 ml

Product characteristics

Filter	yes
Area of the filter	63.6 cm ²
Graduation	yes
Pore size	0.22 µm
Closure type	Screw cap
Connection	Step taper
Example of use	for sterile filtration

Purity & certification

Product category	no medical device IVD
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	12
Piece(s) / case	12
Piece(s) / pallet	144
Height of case	439 mm
Width of case	343 mm
Depth of case	572 mm
Volume of case	0.0861 cbm
Weight of case	4.72 kg
EAN inner box	4038917412148
EAN case	4038917412131

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT