



Product specification

Filtropur V25, Vacuum filtration unit, 250 ml, PES, 0.2 µm



Product information

Product number	83.3940.001
----------------	-------------

Product descriptions

Name	Filtropur V25, Vacuum filtration unit, 250 ml, PES, 0.2 µm
Product number	83.3940.001
Product description	Filtropur V25, Vacuum filtration unit, 250 ml, screw cap, membrane: PES, Ø membraner: 90 mm, pore size: 0.2 µm, for sterile filtration, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/bag

Material & colours

Color of product	transparent
Membrane material	Polyethersulfone (PES)
Closure material	Polypropylene (PP)

Size

Height of product	143.9 mm
Sample volume	250 ml
Filtration volume	250 ml

Product characteristics

Filter	yes
Area of the filter	63.6 cm ²

Graduation	yes
Pore size	0.2 µm
Closure type	Screw cap
Connection	Step taper
Example of use	for sterile filtration

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile, non-pyrogenic/endotoxin-free, non-cytotoxic
Sterilisation	Gamma irradiation
Pyrogen-free	yes
Batched	yes

Packaging

Piece(s) / case	12
Piece(s) / pallet	288
Height of case	381 mm
Width of case	203 mm
Depth of case	521 mm
Volume of case	0.0403 cbm
Weight of case	2.72 kg
EAN inner box	4038917411998
EAN case	4038917411981

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT