




Product specification

**Closure mat, for Deep Well plate 2.2 ml (82.1972.002),
(LxW): 79 x 121 mm, EVA, transparent**



Product information

Product number	95.1991.002
----------------	-------------

Product descriptions

Name	Closure mat, for Deep Well plate 2.2 ml (82.1972.002), (LxW): 79 x 121 mm, EVA, transparent
Product number	95.1991.002
Product description	Closure mat, for Deep Well plate 2.2 ml (82.1972.002), (LxW): 79 x 121 mm, material: EVA, transparent, free from DNA/DNase/RNase, free from pyrogens/endotoxins, 10 piece(s)/bag

Material & colours

Product material	Ethylene-vinyl acetate copolymers (EVA)
Color of product	transparent

Size

Width of product	121 mm
Length of product	79 mm

Product characteristics

Pierceable	no
------------	----

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins
Batched	yes

Packaging

Minimum order quantity	50
Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	10000
Height of case	200 mm
Width of case	280 mm
Depth of case	380 mm
Volume of case	0.0213 cbm
Weight of case	6.4 kg
EAN inner box	4038917135030
EAN case	4038917042031

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT