




Product specification

Closure mat, for Deep Well MegaBlock® 0.5 ml

(82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002), (Lx



Product information

Product number	95.1990.002
----------------	-------------

Product descriptions

Name	Closure mat, for Deep Well MegaBlock® 0.5 ml (82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002), (LxW): 79 x 121 mm, TPE, transparent
Product number	95.1990.002
Product description	Closure mat, for Deep Well MegaBlock® 0.5 ml (82.1969.002) & 1.2 ml (82.1970.002 & 82.1971.002), (LxW): 79 x 121 mm, material: TPE, transparent, pierceable, free from DNA/DNase/RNase, free from pyrogens/endotoxins, 10 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	121 mm
Length of product	79 mm

Product characteristics

Pierceable	yes
------------	-----

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	free from DNA/DNase/RNase, free from pyrogens/endotoxins
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	5500
Height of inner box	65 mm
Width of inner box	205 mm
Depth of inner box	250 mm
Height of case	200 mm
Width of case	280 mm
Depth of case	380 mm
Volume of case	0.0213 cbm
Weight of case	4.1 kg
EAN inner box	4038917135023
EAN case	4038917042048

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT