



## Product specification

**Micro test plate, Terasaki, 60 well, slip-on lid, base shape: conical, PS, transparent**



## Product information

|                |             |
|----------------|-------------|
| Product number | 83.9923.972 |
|----------------|-------------|

## Product descriptions

|                     |   |
|---------------------|---|
| Name                | Micro test plate, Terasaki, 60 well, slip-on lid, base shape: conical, PS, transparent  |
| Product number      | 83.9923.972   |
| Product description | Micro test plate, Terasaki, 60 well, slip-on lid, base shape: conical, lid assembled, material: PS, transparent, sterile, 10 piece(s)/bag |

## Material & colours

|                  |                  |
|------------------|------------------|
| Product material | Polystyrene (PS) |
| Color of product | transparent      |

## Size

|                   |         |
|-------------------|---------|
| Width of product  | 57.8 mm |
| Length of product | 85.1 mm |
| Height of product | 9.7 mm  |
| Max. volume       | 10 µl   |

|                           |             |
|---------------------------|-------------|
| Number of wells           | 60          |
| Base shape                | conical     |
| Vessel type               | Plate       |
| Type of surface treatment | TC-treated  |
| Closure type              | Slip-on lid |

## Purity & certification

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Purity standard  | sterile                 |
| Sterilisation    | Gamma irradiation       |
| Batched          | yes                     |

## Packaging

|                        |               |
|------------------------|---------------|
| Minimum order quantity | 650           |
| Piece(s) / inner box   | 10            |
| Piece(s) / case        | 650           |
| Piece(s) / pallet      | 16250         |
| Height of case         | 264 mm        |
| Width of case          | 324 mm        |
| Depth of case          | 608 mm        |
| Volume of case         | 0.052 cbm     |
| Weight of case         | 12.1998 kg    |
| EAN inner box          | 4038917154710 |
| EAN case               | 4038917354363 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT