



Product specification

Micro test plate, Terasaki, 60 well, slip-on lid, base shape: conical, PS, transparent



Product information

Product number	83.9923.972
----------------	-------------

Product descriptions

Name	Micro test plate, Terasaki, 60 well, slip-on lid, base shape: conical, PS, transparent
Product number	83.9923.972
Product description	Micro test plate, Terasaki, 60 well, slip-on lid, base shape: conical, lid assembled, material: PS, transparent, sterile, 10 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	57.8 mm
Length of product	85.1 mm
Height of product	9.7 mm
Max. volume	10 µl

Product characteristics

Number of wells	60
-----------------	----

Base shape	conical
Vessel type	Plate
Type of surface treatment	TC-treated
Closure type	Slip-on lid

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Gamma irradiation
Batched	yes

Packaging

Piece(s) / inner box	10
Piece(s) / case	650
Piece(s) / pallet	16250
Height of case	264 mm
Width of case	324 mm
Depth of case	608 mm
Volume of case	0.052 cbm
Weight of case	12.1998 kg
EAN inner box	4038917154710
EAN case	4038917354363

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT