



Product specification

**ELISA plate, 96 well, base shape: conical, PS, transparent,
Medium Binding**



Product information

Product number	82.1583.100
----------------	-------------

Product descriptions

Name	ELISA plate, 96 well, base shape: conical, PS, transparent, Medium Binding
Product number	82.1583.100
Product description	ELISA plate, 96 well, base shape: conical, without lid, material: PS, transparent, Medium Binding, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	14.35 mm
Max. volume	290 µl

Product characteristics

Number of wells	96
Base shape	conical
Label/Print	With print
Vessel type	Plate

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order quantity	50
Piece(s) / inner box	25
Piece(s) / case	50
Piece(s) / pallet	2800
Height of case	365 mm
Width of case	232 mm
Depth of case	279 mm
Volume of case	0.0236 cbm
Weight of case	2.64 kg
EAN inner box	4038917356169
EAN case	4038917356152

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT