




## Product specification

### Adhesive film, material: acetate, transparent



## Product information

Product number	82.1586
----------------	---------

## Product descriptions

Name	Adhesive film, Material: Acetate, transparent
Product number	82.1586
Product description	Adhesive film, for ELISA and microtest plates, material: acetate, transparent, self-adhesive, 100 piece(s)/bag

## Material & colours

Color of product	transparent
------------------	-------------

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Batched	yes

## Packaging

Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	30000
Height of inner box	30 mm
Width of inner box	85 mm

Depth of inner box	135 mm
Height of case	225 mm
Width of case	239 mm
Depth of case	95 mm
Volume of case	0.0051 cbm
Weight of case	3 kg
EAN inner box	4038917087148
EAN case	4038917005609

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT