



Product specification

**Lid, for ELISA and microtest plates in 96-well format, PS,
with condensation rings**



Product information

Product number	82.1584
----------------	---------

Product descriptions

Name	Lid, for ELISA and microtest plates in 96-well format, PS, with condensation rings
Product number	82.1584
Product description	Lid, for ELISA and microtest plates in 96-well format, material: PS, with condensation rings, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

Product characteristics

Design	with condensation rings
--------	-------------------------

Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
Batched	yes

Packaging

Piece(s) / inner box	25
Piece(s) / case	100
Piece(s) / pallet	11200
Height of case	235 mm
Width of case	187 mm
Depth of case	282 mm
Volume of case	0.0124 cbm
Weight of case	1.86 kg
EAN inner box	4038917087124
EAN case	4038917005593

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT