



Product specification

**Micro test plate, 96 well, base shape: conical, PS,
transparent**



Product information

Product number	82.1583
----------------	---------

Product descriptions

Name	Micro test plate, 96 Well, Base shape: conical, PS, transparent
Product number	82.1583
Product description	Micro test plate, 96 well, base shape: conical, without lid, material: PS, transparent, free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic, 25 piece(s)/bag

Material & colours

Color of product	transparent
------------------	-------------

Size

Width of product	85.2 mm
Length of product	127.8 mm
Height of product	14.35 mm
Max. volume	290 µl

Product characteristics

Number of wells	96
-----------------	----

Base shape	conical
------------	---------

Vessel type	Plate
-------------	-------

Purity & certification

Product category	not regulated
------------------	---------------

CE certified	no CE mark
--------------	------------

Purity standard	free from DNA/DNase/RNase, free from pyrogens and endotoxins, non-cytotoxic
-----------------	---

Batched	yes
---------	-----

Packaging

Piece(s) / inner box	25
----------------------	----

Piece(s) / case	100
-----------------	-----

Piece(s) / pallet	3200
-------------------	------

Height of case	365 mm
----------------	--------

Width of case	275 mm
---------------	--------

Depth of case	359 mm
---------------	--------

Volume of case	0.036 cbm
----------------	-----------

Weight of case	4.54 kg
----------------	---------

EAN inner box	4038917087117
---------------	---------------

EAN case	4038917355476
----------	---------------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT