




Product specification

Petri dish, 100 x 100 x 20 mm, transparent, without ventilation cams



Product information

Product number	82.9923.422
----------------	-------------

Product descriptions

Name	Petri dish, 100 x 100 x 20 mm, transparent, without ventilation cams
Product number	82.9923.422
Product description	Petri dish, (LxWxH): 100 x 100 x 20 mm, material: PS, transparent, without ventilation cams, sterile, 4 piece(s)/bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Width of product	100 mm
Length of product	100 mm
Height of product	20 mm
Max. volume	175 ml

Product characteristics

Ventilation cams	no
------------------	----

Minimum temperature	0 °C
Maximum temperature	80 °C

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Latex-free	yes
Batched	yes
Safety data sheet	no

Packaging

Minimum order quantity	160
Type of inner packaging	Bag
Piece(s) / case	160
Piece(s) / pallet	4800
Height of case	264 mm
Width of case	324 mm
Depth of case	608 mm
Volume of case	0.052 cbm
Weight of case	6.3501 kg
EAN inner box	4038917154178
EAN case	4038917027854

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT