



## Product specification

### Petri dish, 60 x 15 mm, transparent, with ventilation cams



## Product information

Product number	82.1194.500
----------------	-------------

## Product descriptions

Name	Petri dish, 60 x 15 mm, transparent, with ventilation cams
------	--

Product number	82.1194.500
----------------	-------------

Product description	Petri dish, (ØxH): 60 x 15 mm, material: PS, transparent, with ventilation cams, sterile, 20 piece(s)/tube bag
---------------------	--

## Material & colours

Product material	Polystyrene (PS)
------------------	------------------

Color of product	transparent
------------------	-------------

## Size

Diameter	60 mm
----------	-------

Height of product	15 mm
-------------------	-------

Max. volume	30 ml
-------------	-------

## Product characteristics

Ventilation cams	yes
------------------	-----

Minimum temperature	0 °C
---------------------	------

Maximum temperature	80 °C
---------------------	-------

## Purity & certification

Product category	no medical device   IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Latex-free	yes
Batched	yes
Safety data sheet	no

## Packaging

Minimum order quantity	500
Type of inner packaging	Tube bag
Piece(s) / inner box	20
Piece(s) / case	500
Piece(s) / pallet	12000
Height of case	365 mm
Width of case	317 mm
Depth of case	317 mm
Volume of case	0.0367 cbm
Weight of case	4.49 kg
EAN inner box	4038917347891
EAN case	4038917347884

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT