




Product specification

Petri dish, 150 x 20 mm, transparent, with ventilation cams



Product information

Product number	82.1184.500
----------------	-------------

Product descriptions

Name	Petri dish, 150 x 20 mm, transparent, with ventilation cams
Product number	82.1184.500
Product description	Petri dish, (ØxH): 150 x 20 mm, material: PS, transparent, with ventilation cams, sterile, 10 piece(s)/tube bag

Material & colours

Product material	Polystyrene (PS)
Color of product	transparent

Size

Diameter	150 mm
Height of product	20 mm
Max. volume	275 ml

Product characteristics

Ventilation cams	yes
Minimum temperature	0 °C

Maximum temperature	80 °C
---------------------	-------

Purity & certification

Product category	no medical device IVD
Purity standard	sterile
Sterilisation	Gamma irradiation
Latex-free	yes
Batched	yes
Safety data sheet	no

Packaging

Minimum order quantity	100
Type of inner packaging	Tube bag
Piece(s) / inner box	10
Piece(s) / case	100
Piece(s) / pallet	1800
Height of case	540 mm
Width of case	324 mm
Depth of case	438 mm
Volume of case	0.0766 cbm
Weight of case	6.81 kg
EAN inner box	4038917347860
EAN case	4038917347853

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT