



## Product specification

### Petri dish, 92 x 16 mm, transparent, without ventilation cams



## Product information

Product number	82.1472.001
----------------	-------------

## Product descriptions

Name	Petri dish, 92 x 16 mm, transparent, without ventilation cams
Product number	82.1472.001
Product description	Petri dish, (ØxH): 92 x 16 mm, material: PS, transparent, without ventilation cams, sterile, 20 piece(s)/tube bag

## Material & colours

Color of product	transparent
------------------	-------------

## Size

Diameter	92 mm
Height of product	16 mm
Max. volume	80 ml

## Product characteristics

Ventilation cams	no
Minimum temperature	0 °C
Maximum temperature	80 °C

## Purity & certification

Product category	not regulated
CE certified	no CE mark
Purity standard	sterile
Sterilisation	Gamma irradiation
Latex-free	yes
Batched	yes
Safety data sheet	no

## Packaging

Type of inner packaging	Tube bag
Piece(s) / inner box	480
Piece(s) / case	480
Piece(s) / pallet	9600
Height of case	350 mm
Width of case	387 mm
Depth of case	584 mm
Volume of case	0.0791 cbm
Weight of case	9.1901 kg
EAN inner box	4038917129312
EAN case	4038917022118

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT