



Product specification

Forensic swab, XL, round, in the tube with ventilation membrane, ISO 18385, 125 mm, viscose



Product information

Product number	80.634
----------------	--------

Product descriptions

Name	Forensic swab, XL, round, in the tube with ventilation membrane, ISO 18385, 125 mm, viscose
------	---

Product number	80.634
----------------	--------

Product description	Forensic swab, XL, round, in the tube with ventilation membrane, with safety label, ETO treated, ISO 18385, passive drying, stem length: 125 mm, stem material: PS, swab tip round, swab material: viscose, tube: PP, 50 piece(s)/case
---------------------	--

Material & colours

Color of product	white
------------------	-------

Tube material	Polypropylene (PP)
---------------	--------------------

Stem material	Polystyrene (PS)
---------------	------------------

Swab material	Viscose
---------------	---------

Closure material	Polypropylene (PP)
------------------	--------------------

Color of cap	white
--------------	-------

Length excluding cap	150 mm
Length of product	175 mm
Length of the stem	125 mm

Product characteristics

Design	in the tube with ventilation membrane
Label/Print	With safety label
Cap	Stopper with cotton bud holder
Closure type	Push cap

Purity & certification

Satisfies the requirement	ISO 18385
Product category	not regulated
CE certified	no CE mark
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Type of inner packaging	Box
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	33000
Height of inner box	98 mm
Width of inner box	97 mm
Depth of inner box	181 mm
Height of case	115 mm
Width of case	373 mm
Depth of case	498 mm
Volume of case	0.0214 cbm
Weight of case	4.41 kg
EAN inner box	4038917276702
EAN case	4038917276696

This document was prepared by EDP support
and is valid without signature.



SARSTEDT