



## Product specification

**Forensic swab, round, in the tube with ventilation membrane, ISO 18385, 85 mm, viscose**



## Product information

Product number	80.630
----------------	--------

## Product descriptions

Name	Forensic swab, round, in the tube with ventilation membrane, ISO 18385, 85 mm, viscose
Product number	80.630
Product description	Forensic swab, round, in the tube with ventilation membrane, with safety label, ETO treated, ISO 18385, passive drying, stem length: 85 mm, stem material: PS, swab tip round, swab material: viscose, tube: PP, 50 piece(s)/case

## Material & colours

Color of product	white
Tube material	Polypropylene (PP)
Stem material	Polystyrene (PS)
Swab material	Viscose
Closure material	Polypropylene (PP)
Color of cap	white

Length excluding cap	105 mm
Length of product	175 mm
Length of the stem	85 mm

## Product characteristics

Design	in the tube with ventilation membrane
Label/Print	With safety label
Cap	Standard cap
Closure type	Slip-on lid

## Purity & certification

Satisfies the requirement	ISO 18385
Product category	no medical device   IVD
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Batched	yes

## Packaging

Minimum order quantity	500
Type of inner packaging	Box
Piece(s) / inner box	50
Piece(s) / case	500
Piece(s) / pallet	33000
Height of inner box	98 mm
Width of inner box	97 mm
Depth of inner box	181 mm
Height of case	115 mm
Width of case	373 mm
Depth of case	498 mm
Volume of case	0.0214 cbm
Weight of case	4.52 kg
EAN inner box	4038917275026
EAN case	4038917275019

This document was prepared by EDP support  
and is valid without signature.



# SARSTEDT