




Product specification

Forensic swab, round, ISO 18385, 150 mm, cotton



Product information

Product number	80.626
----------------	--------

Product descriptions

Name	Forensic swab, round, ISO 18385, 150 mm, cotton
Product number	80.626
Product description	Forensic swab, round, ETO treated, ISO 18385, stem length: 150 mm, stem material: wood, swab tip round, swab material: cotton, 1 piece(s)/blister

Material & colours

Color of product	multi-coloured
Stem material	Wood
Swab material	Cotton

Size

Diameter	2.2 mm
Width of product	2.2 mm
Length of product	150 mm
Length of the stem	150 mm

Purity & certification

Satisfies the requirement	ISO 18385
---------------------------	-----------

Product category	not regulated
CE certified	no CE mark
Purity standard	ETO treated
Sterilisation	Gas sterilisation
Batched	yes

Packaging

Type of inner packaging	Bag
Piece(s) / inner box	100
Piece(s) / case	1000
Piece(s) / pallet	68000
Height of inner box	75 mm
Width of inner box	105 mm
Depth of inner box	200 mm
Height of case	135 mm
Width of case	398 mm
Depth of case	428 mm
Volume of case	0.023 cbm
Weight of case	1.95 kg
EAN inner box	4038917279345
EAN case	4038917279338

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT