



Product specification

Swabs, no medium, stem material: PS



Product information

Product number	80.625
----------------	--------

Product descriptions

Name	Swabs, no medium, Stem material: PS
Product number	80.625
Product description	Swabs, in the tube, no medium, (LxØ): 108 x 16 mm, stem material: PS, stem length: 83 mm, swab material: viscose, sterile, 500 piece(s)/bag

Material & colours

Product material	Polypropylene (PP)
Color of product	transparent
Stem material	Polystyrene (PS)
Swab material	Viscose
Closure material	Polyethylene (PE)
Color of cap	transparent

Size

Diameter	16 mm
Length excluding cap	101 mm
Length including cap	108 mm
Length of product	108 mm
Length including cap and plunger	97 mm

Length of the stem

83 mm

Product characteristics

Medium	no medium
Design	in the tube
Label/Print	With paper label
Closure type	Screw cap

Purity & certification

Product category	diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order quantity	500
Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	22000
Height of case	155 mm
Width of case	348 mm
Depth of case	568 mm
Volume of case	0.0306 cbm
Weight of case	3.73 kg
EAN inner box	4038917083003
EAN case	4038917001984

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT