




Product specification

Faeces tube, with blade, screw cap, (LxØ): 55 x 44 mm, white



Product information

Product number	80.9924.027
----------------	-------------

Product descriptions

Name	Faeces tube, with Blade, Screw cap, (LxØ): 55 x 44 mm, white
Product number	80.9924.027
Product description	Faeces tube, with blade, white, screw cap, cap: brown, (LxØ) without cap: 55 x 44 mm, with paper label, 250 piece(s)/bag

Material & colours

Color of product	white
Tube material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Color of cap	brown

Size

Diameter	44 mm
Length excluding cap	55 mm
Length of product	57 mm

Product characteristics

Design	Blade
Label/Print	With paper label
Color of print/label	black
Closure type	Screw cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Batched	yes

Packaging

Piece(s) / inner box	250
Piece(s) / case	500
Piece(s) / pallet	8000
Height of case	373 mm
Width of case	408 mm
Depth of case	458 mm
Volume of case	0.0697 cbm
Weight of case	7.8 kg
EAN inner box	4038917154819
EAN case	4038917039635

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT