



Product specification

Faeces tube, with spoon, push cap, (LxØ): 75 x 23.5 mm, transparent, sterile



Product information

Product number	80.621
----------------	--------

Product descriptions

Name	Faeces tube, with Spoon, Push cap, (LxØ): 75 x 23.5 mm, transparent, sterile
Product number	80.621
Product description	Faeces tube, with spoon, transparent, push cap, cap: brown, (LxØ) without cap: 75 x 23.5 mm, sterile, 50 piece(s)/bag

Material & colours

Color of product	transparent
Tube material	Polystyrene (PS)
Closure material	Low Density Polyethylene (LD-PE)
Color of cap	brown

Size

Diameter	23.5 mm
Length excluding cap	75 mm
Length of product	87 mm

Product characteristics

Design	Spoon
Closure type	Push cap

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Piece(s) / inner box	50
Piece(s) / case	250
Piece(s) / pallet	13000
Height of inner box	60 mm
Width of inner box	305 mm
Depth of inner box	350 mm
Height of case	130 mm
Width of case	380 mm
Depth of case	548 mm
Volume of case	0.0271 cbm
Weight of case	2.73 kg
EAN inner box	4038917082969
EAN case	4038917001946

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT