



Product specification

SAHARA-TSC 115 V



Product information

Product number	97.8710.602
----------------	-------------

Product descriptions

Name	SAHARA-TSC 115 V
------	------------------

Product number	97.8710.602
----------------	-------------

Product description	SAHARA-TSC 115 V, for thawing and warming blood products in plastic bags before transfusion, suitable for cryopreserved stem cells, incl. adaptation compress and warming tray
---------------------	--

Size

Width of product	32 cm
------------------	-------

Length of product	49.3 cm
-------------------	---------

Height of product	32.5 cm
-------------------	---------

Product characteristics

Power consumption	655 W
-------------------	-------

Supply voltage	115 V
----------------	-------

Purity & certification

Product category	Medical device equipment, CE
------------------	------------------------------

CE certified	CE - manufacturer's self-declaration
--------------	--------------------------------------

Batched	no
---------	----

Packaging

Volume of case	0.1639 cbm
Weight of case	23.75 kg
EAN inner box	4038917968577
EAN case	4038917968584

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



SARSTEDT