



# Product specification

## Adaptation compress TSC



### Product information

Product number	79.8710.610
----------------	-------------

### Product descriptions

Name	Adaptation compress TSC
Product number	79.8710.610
Product description	Adaptation compress TSC, Replacement gel compress for thawing out cryopreserved blood stem cell preparations in the SAHARA TSC

### Size

Width of product	16 cm
Length of product	23.5 cm

### Purity & certification

Product category	no medical device   IVD
Batched	yes

### Packaging

Piece(s) / pallet	960
Height of case	140 mm
Width of case	220 mm
Depth of case	320 mm
Volume of case	0.0099 cbm

Weight of case	2.6 kg
EAN inner box	4038917978545
EAN case	4038917978552

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

This document was prepared by EDP support and is valid without signature.



# SARSTEDT